

Deluxe Comfort self propelled wheelchair - Standard



Product Description

A modern looking self propelled wheelchair with many features. Has attendant brakes and solid wheels to prevent punctures. The large wheels allow the user to move themselves and lighter if pushed by a carer. The wheels will also accommodate a smoother ride and is better suited for rougher terrain than an attendant wheelchair. The arm rests adjust to suit comfort of the user and for a better leverage position for help standing up from the wheelchair. Has quick and easy release rear wheels to fold more compact and lighter when getting into a car. Comes with Velcro fashioning back and seat cushions. Comes in seat width choice of 18" and 20". Has swing away footrests that are height adjustable. Colour Silver.

Product Code: WC0303

£249.99

£299.99 inc VAT

You may be able to claim VAT relief.

Specification

Weight:	14kg	User Weight Limit:	115kg
Height:	87cm	Width:	87cm
Depth:	103cm	Seat Height:	50cm
Seat Width:	44cm	Seat Depth:	42cm
Handle Height:	87cm	Colour:	Silver
Warranty:	1 Year		

All Prices Correct on 02/Jul/2022.

Sales Support: 01743 250 821
 General Telephone: 01743 250 820
 Email: info@ilp-ltd.co.uk
 Website: www.ilp-ltd.co.uk

The Lantern, Meadow Farm Drive, Shrewsbury, Shropshire, SY1 4NG.

