

Handybar Car Transfer Handle



Product Code: DL0506

£34.99

£41.99 inc VAT

You may be able to claim VAT relief.

Product Description

The Handy Bar is a great tool designed for car use. The top part of the bar is designed to slot into a car doors locking mechanism, once in place this then secures the handle part into a convenient grab rail. The bar can be placed in 2 different positions and fits most car doors. The Handy Bar also features an emergency seat belt cutter, the blade is integral to the bar and is not possible to reach with fingertips. The top of the bar is hardened steel and can be used as a hammer to break the windscreen in an emergency.

Specification

Weight:	0kg	User Weight Limit:	kg
Height:	23cm	Width:	cm
Depth:	cm	Seat Height:	cm
Seat Width:	cm	Seat Depth:	cm
Handle Height:	cm	Colour:	Red
Warranty:	None		

All Prices Correct on 16/Aug/2022.

Sales Support: 01743 250 821
 General Telephone: 01743 250 820
 Email: info@ilp-ltd.co.uk
 Website: www.ilp-ltd.co.uk

The Lantern, Meadow Farm Drive, Shrewsbury, Shropshire, SY1 4NG.

