

Topro Neuro - Medium Grey



Product Code: MB0662

£329.99

£395.99 inc VAT

You may be able to claim VAT relief.

Product Description

The Neuro has a reverse braking system for extra safety measures. The brakes of the walker remain on until the user squeezes the brake handles to release them. This ensures the person with advanced mobility problems can walk slowly and securely safe with the knowledge the brakes are easily applied by letting go of the brake handle. Furthermore the seat is safe to use at any time without worrying about locking brakes beforehand. The Neuro folds into an x shape and the handles can be adjusted down for storage into the car. The Neuro comes in 2 different sizes, small and medium. Small has a seat height of 54cm with handles adjustable between 71-86cm. The Medium has a seat height of 62cm with handles adjustable between 78-100cm. The Neuro also comes with a free back support.

Specification

Weight:	7kg	User Weight Limit:	150kg
Height:	100cm	Width:	60cm
Depth:	65cm	Seat Height:	62cm
Seat Width:	47cm	Seat Depth:	cm
Handle Height:	100cm	Colour:	Grey
Warranty:	7 Years		

All Prices Correct on 27/May/2022.

Sales Support:

01743 250 821

General Telephone:

01743 250 820

Email:

info@ilp-ltd.co.uk

Website:

www.ilp-ltd.co.uk

The Lantern, Meadow Farm Drive, Shrewsbury, Shropshire, SY1 4NG.

